



SILABE

Samples request form
mail to : samples@silabe.com

DATE:

Company

Contact:

Phone:

Address:

City, State, Zip:

Email:

Macaca mulatta - Macaque rhesus - Rhesus macaque

		Volume (ml)	Qty	Specific gender	Temperature
Whole blood	□ LH	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
	□ K2-EDTA	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
	□ K3-EDTA	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
□ NH	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
	□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
□ Na Citrate	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
	□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
□ CTAD	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
	□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
CSF	□ CSF	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen

		Volume (ml)	Qty	Specific gender	Temperature
Plasma	□ LH	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
	□ K2-EDTA	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
	□ K3-EDTA	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
□ NH	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
	□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
□ Na Citrate	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
	□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen	
Serum	□ Serum	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
Bone marrow	□ Bone marrow	□ Ind.tube		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen
		□ Pool		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen

		Qty of CPT per animal	Animal number	Specific gender	Temperature
PBMC	□ PBMC	□ CPT tube/Animal		□ M ____ □ F ____ □ Indifferent	□ -80°C

CPT tube	□ CPT tube	centrifugated □ Yes □ No		□ M ____ □ F ____ □ Indifferent	□ Room T° □ Chilled T° □ Frozen

			Temperature
Others	□ Others non sterile		
	□ Feces		□ Chilled T° □ Frozen
	□ Urine		□ Chilled T° □ Frozen

Comments